

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

In re application of:

YU *et al.*

Appl. No. 09/314,889

Filed: May 19, 1999

For: **Death Domain Containing
Receptors**

Art Unit: 1646

Examiner: Ulm, J.

Atty. Docket: 1488.0310006/EKS/SGW

Declaration for Deposited Biological Materials

Commissioner for Patents
Washington, D.C. 20231


Sir:

The undersigned attorney of record states as follows:

1. A plasmid containing human cDNA encoding DR3-V1 protein was deposited under the terms of the Budapest Treaty on October 10, 1996. The deposit was made at the American Type Culture Collection (ATCC), 10801 University Boulevard, Manassas, Virginia 20110-2209 (current address), and given accession number 97757. A copy of the ATCC deposit receipt reciting the term of the deposit, as well as the statement of viability, is attached hereto.

2. Assurance is hereby given that all restrictions on the availability to the public of the deposited plasmid referred to above will be irrevocably removed upon the granting of a patent issuing from the patent application captioned above.

DECEMBER 4, 2000
Date


Jonathan L. Klein
Attorney for Applicant
Registration No. 41,119

Human Genome Sciences, Inc.
9410 Key West Avenue
Rockville, MD 20850
Telephone: (301) 251-6015

031-6.DRC.DEPOSIT.WPD



American Type Culture Collection

12301 Parklawn Drive • Rockville, MD 20852 USA • Telephone: (301) 231-5528 Telex: 898-055 ATCCNORTH • FAX: 301-770-2587

BUDAPEST TREATY ON THE INTERNATIONAL RECOGNITION OF THE DEPOSIT OF MICROORGANISMS FOR THE PURPOSES OF PATENT PROCEDURE

INTERNATIONAL FORM

RECEIPT IN THE CASE OF AN ORIGINAL DEPOSIT ISSUED PURSUANT TO RULE 7.3
AND VIABILITY STATEMENT ISSUED PURSUANT TO RULE 10.2

MAR 13 1996
HGS PATENT DEPT.

To: (Name and Address of Depositor or Attorney)

Human Genome Sciences, Inc.
Attn: Robert H. Benson
9410 Key West Avenue
Rockville, MD 20850

Deposited on Behalf of: Human Genome Sciences, Inc.

Identification Reference by Depositor:

ATCC Designation

DNA Plasmid, 231556

97456 ~ .

The deposits were accompanied by: a scientific description a proposed taxonomic description indicated above.

The deposits were received March 1, 1996 by this International Depository Authority and have been accepted.

AT YOUR REQUEST:

☒ We will inform you of requests for the strains for 30 years.

The strains will be made available if a patent office signatory to the Budapest Treaty certifies one's right to receive, or if a U.S. Patent is issued citing the strains, and ATCC is instructed by the United States Patent & Trademark Office or the depositor to release said strains.

If the cultures should die or be destroyed during the effective term of the deposit, it shall be your responsibility to replace them with living cultures of the same.

The strains will be maintained for a period of at least 30 years from date of deposit, or five years after the most recent request for a sample, whichever is longer. The United States and many other countries are signatory to the Budapest Treaty.

The viability of the cultures cited above was tested March 7, 1996. On that date, the cultures were viable.

International Depository Authority: American Type Culture Collection, Rockville, Md. 20852 USA

Signature of person having authority to represent ATCC:

Barbara M. Hailey
Barbara M. Hailey, Administrator, Patent Depository

Date: March 8, 1996



American Type Culture Collection

12301 Parklawn Drive • Rockville, MD 20852 USA • Telephone: (301) 231-5520 Telex: 908-768 ATCCROVE • FAX: 301-816-4344

BUDAPEST TREATY ON THE INTERNATIONAL RECOGNITION OF THE DEPOSIT OF MICROORGANISMS FOR THE PURPOSES OF PATENT PROCEDURE

INTERNATIONAL FORM

RECEIPT IN THE CASE OF AN ORIGINAL DEPOSIT ISSUED PURSUANT TO RULE 7.3
AND VIABILITY STATEMENT ISSUED PURSUANT TO RULE 10.2

To: (Name and Address of Depositor or Attorney)

Human Genome Sciences, Inc.
Attn: Robert H. Benson
9410 Key West Avenue
Rockville, MD 20850

Deposited on Behalf of: Human Genome Sciences, Inc.

Identification Reference by Depositor:

ATCC Designation

DNA Plasmid, 1198782 (Docket PF267P1.SKB)

97757

The deposits were accompanied by: ☐ a scientific description ☐ a proposed taxonomic description indicated above.

The deposits were received October 10, 1996 by this International Depository Authority and have been accepted.

AT YOUR REQUEST:

☒ We will inform you of requests for the strains for 30 years.

The strains will be made available if a patent office signatory to the Budapest Treaty certifies one's right to receive, or if a U.S. Patent is issued citing the strains, and ATCC is instructed by the United States Patent & Trademark Office or the depositor to release said strains.

If the cultures should die or be destroyed during the effective term of the deposit, it shall be your responsibility to replace them with living cultures of the same.

The strains will be maintained for a period of at least 30 years from date of deposit, or five years after the most recent request for a sample, whichever is longer. The United States and many other countries are signatory to the Budapest Treaty.

The viability of the cultures cited above was tested October 17, 1996. On that date, the cultures were viable.

International Depository Authority: American Type Culture Collection, Rockville, Md. 20852 USA

Signature of person having authority to represent ATCC:

Barbara M. Hailey, Administrator, Patent Depository

Date: October 19, 1996

CURRENT PROTOCOLS IN MOLECULAR BIOLOGY

VOLUME 1

EDITORIAL BOARD

Frederick M. Ausubel

Massachusetts General Hospital & Harvard Medical School

Roger Brent

Massachusetts General Hospital & Harvard Medical School

Robert E. Kingston

Massachusetts General Hospital & Harvard Medical School

David D. Moore

Massachusetts General Hospital & Harvard Medical School

J.G. Seidman

Harvard Medical School

John A. Smith

University of Alabama

Kevin Struhl

Harvard Medical School

GUEST EDITORS

Lisa M. Albright

DNA Sequencing

Donald M. Coen

Harvard Medical School

Polymerase Chain Reaction

Ajit Varki

University of California San Diego

Glycoproteins

SERIES EDITOR

Virginia Benson Chanda



John Wiley & Sons, Inc.

CORE 13 (S36)

Copyright © 1994–1997 by John Wiley & Sons, Inc.

Copyright © 1987–1994 by Current Protocols

All rights reserved. Published simultaneously in Canada.

Reproduction or translation of any part of this work beyond that permitted by Section 107 or 108 of the 1976 United States Copyright Act without the permission of the copyright owner is unlawful. Requests for permission or further information should be addressed to the Permissions Department, John Wiley & Sons, Inc.

While the authors, editors, and publisher believe that the specification and usage of reagents, equipment, and devices, as set forth in this book, are in accord with current recommendations and practice at the time of publication, they accept no legal responsibility for any errors or omissions, and make no warranty, express or implied, with respect to material contained herein. In view of ongoing research, equipment modifications, changes in governmental regulations, and the constant flow of information relating to the use of experimental reagents, equipment, and devices, the reader is urged to review and evaluate the information provided in the package insert or instructions for each chemical, piece of equipment, reagent, or device for, among other things, any changes in the instructions or indication of usage and for added warnings and precautions. This is particularly important in regard to new or infrequently employed chemicals or experimental reagents.

Library of Congress Cataloging in Publication Data:

Current protocols in molecular biology. 3 vols.

1. Molecular biology—Technique. 2. Molecular biology—Laboratory manuals. I. Ausubel, Frederick M.

QH506.C87 1987 574.8'8'028 87-21033
ISBN 0-471-50338-X

Printed in the United States of America

20 19 18 17 16 15 14 13.

IN SITU HYBRIDIZATION AND IMMUNOHISTOCHEMISTRY

14

INTRODUCTION	14.0.3
14.1 Fixation, Embedding, and Sectioning of Tissues, Embryos, and Single Cells	14.1.1
Basic Protocol: Paraformaldehyde Fixation and Paraffin Wax Embedding of Tissues and Embryos	14.1.1
Alternate Protocol: Fixation of Suspended and Cultured Cells	14.1.3
Support Protocol: Perfusion of Adult Mice	14.1.4
Support Protocol: Sectioning Samples in Wax Blocks	14.1.4
Reagents and Solutions	14.1.6
Commentary	14.1.7
14.2 Cryosectioning	14.2.1
Basic Protocol: Specimen Preparation and Sectioning	14.2.1
Support Protocol: Fixation of Cryosections for In situ Hybridization	14.2.5
Support Protocol: Tissue Fixation and Sucrose Infusion	14.2.6
Commentary	14.2.6
14.3 In situ Hybridization to Cellular RNA	14.3.1
Basic Protocol: Hybridization Using Paraffin Sections and Cells	14.3.1
Alternate Protocol: Hybridization Using Cryosections	14.3.5
Support Protocol: Synthesis of ³⁵ S-Labeled Riboprobes	14.3.7
Support Protocol: Synthesis of ³⁵ S-Labeled Double-Stranded DNA Probes	14.3.8
Reagents and Solutions	14.3.8
Commentary	14.3.11
14.4 Detection of Hybridized Probe	14.4.1
Basic Protocol: Film Autoradiography	14.4.1
Basic Protocol: Emulsion Autoradiography	14.4.1
Support Protocol: Preparation of Diluted Emulsion for Autoradiography	14.4.2
Commentary	14.4.3
14.5 Counterstaining and Mounting of Autoradiographed In situ Hybridization Slides	14.5.1
Basic Protocol: Giemsa Staining	14.5.1
Alternate Protocol: Hematoxylin/Eosin Staining	14.5.2
Alternate Protocol: Toluidine Blue Staining	14.5.3
Alternate Protocol: Hoechst Staining	14.5.3
Reagents and Solutions	14.5.4
Commentary	14.5.5
14.6 Immunohistochemistry	14.6.1
Basic Protocol: Immunofluorescent Labeling of Cells Grown as Monolayers	14.6.1
Alternate Protocol: Immunofluorescent Labeling of Suspension Cells	14.6.2
Basic Protocol: Immunofluorescent Labeling of Tissue Sections	14.6.3
Alternate Protocol: Immunofluorescent Labeling of Tissue Sections Using Coplin Jars	14.6.4
Alternate Protocol: Immunofluorescent Labeling Using Streptavidin-Biotin Conjugates	14.6.5
Alternate Protocol: Immunogold Labeling of Tissue Sections	14.6.6
Alternate Protocol: Immunoperoxidase Labeling of Tissue Sections	14.6.6
Alternate Protocol: Immunofluorescent Double-Labeling of Tissue Sections	14.6.7
Reagents and Solutions	14.6.7
Commentary	14.6.8

continued

14.0.1

SDS electrophoresis buffer, 5×

15.1 g Tris base

72.0 g glycine

5.0 g SDS

H₂O to 1000 ml

Dilute to 1× or 2× for working solution, as appropriate

Do not adjust the pH of the stock solution, as the solution is pH 8.3 when diluted. Store at 0° to 4°C until use (up to 1 month).

SED (standard enzyme diluent)

20 mM Tris·Cl, pH 7.5

500 µg/ml bovine serum albumin (Pentax Fraction V)

10 mM 2-mercaptoethanol

Store up to 1 month at 4°C

Sodium acetate, 3 M

Dissolve 408 g sodium acetate·3H₂O in 800 ml H₂O

Add H₂O to 1 liter

Adjust pH to 4.8 or 5.2 (as desired) with 3 M acetic acid

Sodium acetate buffer, 0.1 M

Solution A: 11.55 ml glacial acetic acid/liter (0.2 M).

Solution B: 27.2 g sodium acetate (NaC₂H₃O₂·3H₂O)/liter (0.2 M).

Referring to Table A.2.2 for desired pH, mix the indicated volumes of solutions A and B, then dilute with H₂O to 100 ml. (See Potassium acetate buffer recipe for further details.)

Sodium phosphate buffer, 0.1 M

Solution A: 27.6 g NaH₂PO₄·H₂O per liter (0.2 M).

Solution B: 53.65 g Na₂HPO₄·7H₂O per liter (0.2 M).

Referring to Table A.2.3 for desired pH, mix the indicated volumes of solutions A and B, then dilute with H₂O to 200 ml. (See Potassium phosphate buffer recipe for further details.)

SSC (sodium chloride/sodium citrate), 20×

3 M NaCl (175 g/liter)

0.3 M Na₃citrate·2H₂O (88 g/liter)

Adjust pH to 7.0 with 1 M HCl

STE buffer

10 mM Tris·Cl, pH 7.5

10 mM NaCl

1 mM EDTA, pH 8.0

TAE (Tris/acetate/EDTA) electrophoresis buffer

50× stock solution:

Working solution, pH ~8.5:

242 g Tris base

57.1 ml glacial acetic acid

37.2 g Na₂EDTA·2H₂O

H₂O to 1 liter

40 mM Tris·acetate

2 mM Na₂EDTA·2H₂O

TBE (Tris/borate/EDTA) electrophoresis buffer

10× stock solution, 1 liter:

108 g Tris base (890 mM)

55 g boric acid (890 mM)

40 ml 0.5 M EDTA, pH 8.0 (20 mM)

CURRENT PROTOCOLS IN MOLECULAR BIOLOGY

VOLUME 1

EDITORIAL BOARD

Frederick M. Ausubel

Massachusetts General Hospital & Harvard Medical School

Roger Brent

Massachusetts General Hospital & Harvard Medical School

Robert E. Kingston

Massachusetts General Hospital & Harvard Medical School

David D. Moore

Massachusetts General Hospital & Harvard Medical School

J.G. Seidman

Harvard Medical School

John A. Smith

University of Alabama

Kevin Struhl

Harvard Medical School

GUEST EDITORS

Lisa M. Albright

DNA Sequencing

Donald M. Coen

Harvard Medical School
Polymerase Chain Reaction

Ajit Varki

University of California San Diego
Glycoproteins

SERIES EDITOR

Virginia Benson Chanda



John Wiley & Sons, Inc.

CORE 13 (S36)

Copyright © 1994–1997 by John Wiley & Sons, Inc.

Copyright © 1987–1994 by Current Protocols

All rights reserved. Published simultaneously in Canada.

Reproduction or translation of any part of this work beyond that permitted by Section 107 or 108 of the 1976 United States Copyright Act without the permission of the copyright owner is unlawful. Requests for permission or further information should be addressed to the Permissions Department, John Wiley & Sons, Inc.

While the authors, editors, and publisher believe that the specification and usage of reagents, equipment, and devices, as set forth in this book, are in accord with current recommendations and practice at the time of publication, they accept no legal responsibility for any errors or omissions, and make no warranty, express or implied, with respect to material contained herein. In view of ongoing research, equipment modifications, changes in governmental regulations, and the constant flow of information relating to the use of experimental reagents, equipment, and devices, the reader is urged to review and evaluate the information provided in the package insert or instructions for each chemical, piece of equipment, reagent, or device for, among other things, any changes in the instructions or indication of usage and for added warnings and precautions. This is particularly important in regard to new or infrequently employed chemicals or experimental reagents.

Library of Congress Cataloging in Publication Data:

Current protocols in molecular biology. 3 vols.

1. Molecular biology—Technique. 2. Molecular biology—Laboratory manuals. I. Ausubel, Frederick M.

QH506.C87 1987 574.8'8'028 87-21033

ISBN 0-471-50338-X

Printed in the United States of America

20 19 18 17 16 15 14 13

IN SITU HYBRIDIZATION AND IMMUNOHISTOCHEMISTRY

14

INTRODUCTION	14.0.3
14.1 Fixation, Embedding, and Sectioning of Tissues, Embryos, and Single Cells	14.1.1
Basic Protocol: Paraformaldehyde Fixation and Paraffin Wax Embedding of Tissues and Embryos	14.1.1
Alternate Protocol: Fixation of Suspended and Cultured Cells	14.1.3
Support Protocol: Perfusion of Adult Mice	14.1.4
Support Protocol: Sectioning Samples in Wax Blocks	14.1.4
Reagents and Solutions	14.1.6
Commentary	14.1.7
14.2 Cryosectioning	14.2.1
Basic Protocol: Specimen Preparation and Sectioning	14.2.1
Support Protocol: Fixation of Cryosections for In situ Hybridization	14.2.5
Support Protocol: Tissue Fixation and Sucrose Infusion	14.2.6
Commentary	14.2.6
14.3 In situ Hybridization to Cellular RNA	14.3.1
Basic Protocol: Hybridization Using Paraffin Sections and Cells	14.3.1
Alternate Protocol: Hybridization Using Cryosections	14.3.5
Support Protocol: Synthesis of ³⁵ S-Labeled Riboprobes	14.3.7
Support Protocol: Synthesis of ³⁵ S-Labeled Double-Stranded DNA Probes	14.3.8
Reagents and Solutions	14.3.8
Commentary	14.3.11
14.4 Detection of Hybridized Probe	14.4.1
Basic Protocol: Film Autoradiography	14.4.1
Basic Protocol: Emulsion Autoradiography	14.4.1
Support Protocol: Preparation of Diluted Emulsion for Autoradiography	14.4.2
Commentary	14.4.3
14.5 Counterstaining and Mounting of Autoradiographed In situ Hybridization Slides	14.5.1
Basic Protocol: Giemsa Staining	14.5.1
Alternate Protocol: Hematoxylin/Eosin Staining	14.5.2
Alternate Protocol: Toluidine Blue Staining	14.5.3
Alternate Protocol: Hoechst Staining	14.5.3
Reagents and Solutions	14.5.4
Commentary	14.5.5
14.6 Immunohistochemistry	14.6.1
Basic Protocol: Immunofluorescent Labeling of Cells Grown as Monolayers	14.6.1
Alternate Protocol: Immunofluorescent Labeling of Suspension Cells	14.6.2
Basic Protocol: Immunofluorescent Labeling of Tissue Sections	14.6.3
Alternate Protocol: Immunofluorescent Labeling of Tissue Sections Using Coplin Jars	14.6.4
Alternate Protocol: Immunofluorescent Labeling Using Streptavidin-Biotin Conjugates	14.6.5
Alternate Protocol: Immunogold Labeling of Tissue Sections	14.6.6
Alternate Protocol: Immunoperoxidase Labeling of Tissue Sections	14.6.6
Alternate Protocol: Immunofluorescent Double-Labeling of Tissue Sections	14.6.7
Reagents and Solutions	14.6.7
Commentary	14.6.8

continued

- 1c. **Harsh treatment:** Pour several hundred milliliters of boiling 0.1% SDS onto the membrane. Cool to room temperature.

If a membrane is to be reprobed, it must not be allowed to dry out between hybridization and stripping. If it becomes dry, the probe may bind to the matrix.

2. Place membrane on a sheet of dry Whatman 3MM filter paper and blot excess liquid with a second sheet. Wrap the membrane in plastic wrap and set up an autoradiograph.

If signal is still seen after autoradiography, rewash using harsher conditions.

3. The membrane can now be rehybridized. Alternatively, it can be dried and stored for later use.

Membranes can be stored dry between Whatman 3MM paper for several months at room temperature. For long-term storage, place the membranes in a desiccator at room temperature or 4°C.

REAGENTS AND SOLUTIONS

Aqueous prehybridization/hybridization (APH) solution

5× SSC (APPENDIX 2)

5× Denhardt solution (APPENDIX 2)

1% (w/v) SDS

Add 100 µg/ml denatured salmon sperm DNA (see below) just before use

Alternatives to Denhardt solution and denatured salmon sperm DNA as blocking agents are listed in Table 2.10.5 (see discussion in critical parameters).

Denatured salmon sperm DNA

Dissolve 10 mg Sigma type III salmon sperm DNA (sodium salt) in 1 ml water. Pass vigorously through a 17-G needle 20 times to shear the DNA. Place in a boiling water bath for 10 min, then chill. Use immediately or store at -20°C in small aliquots. If stored, reheat to 100°C for 5 min and chill on ice immediately before using.

Formamide prehybridization/hybridization (FPH) solution

5× SSC (APPENDIX 2)

5× Denhardt solution (APPENDIX 2)

50% (w/v) formamide

1% (w/v) SDS

Add 100 µg/ml denatured salmon sperm DNA (see above) just before use

Alternatives to Denhardt solution and denatured salmon sperm DNA as blocking agents are listed in Table 2.10.5 (see discussion in critical parameters).

Commercial formamide is usually satisfactory for use. If the liquid has a yellow color, deionize as follows: add 5 g of mixed-bed ion-exchange resin [e.g., Bio-Rad AG 501-X8 or 501-X8(D) resins] per 100 ml formamide, stir at room temperature for 1 hr, and filter through Whatman no. 1 paper.

CAUTION: Formamide is a teratogen. Handle with care.

Labeling buffer

200 mM Tris·Cl, pH 7.5

30 mM MgCl₂

10 mM spermidine

Mild stripping solution

5 mM Tris·Cl, pH 8.0

2 mM EDTA

0.1× Denhardt solution (APPENDIX 2)